

## Items Needed:

- Item 1:** A container that is big enough for handpieces to sit in completely, with corresponding lid.
- Item 2:** A zip lock bag big enough for handpieces to fit in and be securely closed.
- Item 3:** Enough Lubricant to fully submerge handpiece in container.
- Item 4:** A secure place to keep handpiece.
- 

## Lubricate & Soak Handpiece/Air Motor:

**Step 1: Detach Handpiece/Air Motor** from unit tubing.

**Step 2:** Clean Handpiece/Air Motor after use.

**Step 3: Soak Handpiece/Air Motor in Lubricant.** Using a small tupperware container with a lid, pour lubricant into the container submerging the air motor and placing the lid back onto the container. **Soak the Handpiece/Air Motor over night.**

**Step 4: The Next Day,** remove Handpiece/Air Motor from the lubricant and place into a zip lock bag. **Do Not Sterilize the Handpiece/Air Motor Just Yet. Let the lubricant stay in the Handpiece/ Air Motor. Leave as is, until Handpiece/Air Motor is needed.**

**Step 5: The Day Before the Handpiece/Air Motor is needed, Sterilize Handpiece/Air Motor.**

**Step 6:** After use of Handpiece/Air Motor, **Repeat Steps 1 through 5.**



We Recommend:

**Soaking your Handpieces in Lubricant Overnight; at Least Once a Month.**

*This will help Remove Stuck on Debris & Saturate Bearings with Lubricant, for longer life.*

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

---

Helpful Links:

We have lots of extremely useful information on our website:

[www.mcshandpiece.com](http://www.mcshandpiece.com)

For an In-depth Look at Handpiece Maintenance:

[www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

Cause & Effect of Excessive Pressure Applied to the Bur:

[www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur](http://www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur)

Understanding High Speed Bearings:

[www.mcshandpiece.com/understanding-high-speed-bearings](http://www.mcshandpiece.com/understanding-high-speed-bearings)

How to Free Spin Your Turbine:

[www.mcshandpiece.com/free-spin](http://www.mcshandpiece.com/free-spin)

Handpiece Lubricating Procedure:

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

Chuck Mechanism Cleaning Procedure:

[www.mcshandpiece.com/chuck-brush-procedure](http://www.mcshandpiece.com/chuck-brush-procedure)

Specialized Maintenance Procedures:

Highspeed's: [www.mcshandpiece.com/highspeed-maintenance](http://www.mcshandpiece.com/highspeed-maintenance)

Low Speed's: [www.mcshandpiece.com/low-speed-maintenance](http://www.mcshandpiece.com/low-speed-maintenance)

Electric Air Powered Highspeed's: [www.mcshandpiece.com/electric-airpowered-maintenance](http://www.mcshandpiece.com/electric-airpowered-maintenance)

Air Motors: [www.mcshandpiece.com/air-motor-maintenance](http://www.mcshandpiece.com/air-motor-maintenance)

Accessories & Other Purchases:

[www.mcshandpiece.com/products](http://www.mcshandpiece.com/products)

---

For More Information about Maintenance: [www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

T: (604) 575-8544 - TF: (888) 575-8544 - TXT: (604) 373-6747 - [info@mcshandpiece.com](mailto:info@mcshandpiece.com) - [www.mcshandpiece.com](http://www.mcshandpiece.com)

**Office Hours:** Monday – Friday; 8am – 2pm (PST) - **Office Closures:** Please Visit Our Website for Dates!

Page 2 of 2