

## Clean & Chuck Mechanism.

**Step 1: Lubricate Chuck Mechanism with 1 Drop of Lubricant.**

*(Using Health Canada/FDA Approved Mineral Oil Lubricant)*



**Step 2: Insert Chuck Brush into Chuck Mechanism.**



**Step 3: Turn Chuck Brush in a Circular Motion, Back & Forth for 5 seconds.**

**Step 4: Lubricate Chuck Mechanism with 1 Drop of Lubricant, Repeating Step 1.**

**Step 5: Run Handpiece (With Bur or Shipping Post) for 5 seconds to expel debris.**

**Chuck Brush Is Reusable!  
Please Ensure Suitable Cleaning of Brush after Every Use.**

[www.mcshandpiece.com/accessories](http://www.mcshandpiece.com/accessories)

## **Lubricate & Clean the Handpiece.**

**Step 1:** De-attach **Handpiece** from the **Air Motor**.

**Step 2:** Wipe off parts with **Alcohol**.

**Step 3:** Hold the **Handpiece Upside Down**. The **Opening** for the **Air Motor** should be **Facing Upward**. With an **Aerosol Lubricant Spray** (Shown in Photo to Left) & **Fill the Cavity** with **Lubricant**.



**Step 4:** While Still Holding the Handpiece Upside Down, Attach the Motor & the unit tubing to the handpiece. Run for 30 Seconds.

**Step 5:** After **Autoclaving the Handpiece**, Allow to Cool Gradually for **30 Minutes** to **Room Temperature**.

**Step 6:** Repeat Steps 2 - 4.

## We Recommend:

**Soaking your Handpieces in Lubricant Overnight; at Least Once a Month.**

*This will help Remove Stuck on Debris & Saturate Bearings with Lubricant, for longer life.*

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

---

## Helpful Links:

We have lots of extremely useful information on our website:

[www.mcshandpiece.com](http://www.mcshandpiece.com)

For an In-depth Look at Handpiece Maintenance:

[www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

Cause & Effect of Excessive Pressure Applied to the Bur:

[www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur](http://www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur)

Understanding High Speed Bearings:

[www.mcshandpiece.com/understanding-high-speed-bearings](http://www.mcshandpiece.com/understanding-high-speed-bearings)

How to Free Spin Your Turbine:

[www.mcshandpiece.com/free-spin](http://www.mcshandpiece.com/free-spin)

Handpiece Lubricating Procedure:

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

Chuck Mechanism Cleaning Procedure:

[www.mcshandpiece.com/chuck-brush-procedure](http://www.mcshandpiece.com/chuck-brush-procedure)

Specialized Maintenance Procedures:

Highspeed's: [www.mcshandpiece.com/highspeed-maintenance](http://www.mcshandpiece.com/highspeed-maintenance)

Low Speed's: [www.mcshandpiece.com/low-speed-maintenance](http://www.mcshandpiece.com/low-speed-maintenance)

Electric Air Powered Highspeed's: [www.mcshandpiece.com/electric-airpowered-maintenance](http://www.mcshandpiece.com/electric-airpowered-maintenance)

Air Motors: [www.mcshandpiece.com/air-motor-maintenance](http://www.mcshandpiece.com/air-motor-maintenance)

Accessories & Other Purchases:

[www.mcshandpiece.com/products](http://www.mcshandpiece.com/products)

---

For More Information about Maintenance: [www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

T: (604) 575-8544 - TF: (888) 575-8544 - TXT: (604) 373-6747 - [info@mcshandpiece.com](mailto:info@mcshandpiece.com) - [www.mcshandpiece.com](http://www.mcshandpiece.com)

**Office Hours:** Monday – Friday; 8am – 2pm (PST) - **Office Closures:** Please Visit Our Website for Dates!

Page 3 of 3