

## Clean the Chuck Mechanism.

**Step 1: Lubricate Chuck Mechanism with 1 Drop of Lubricant.**

*(Using Health Canada/FDA Approved Mineral Oil Lubricant)*



**Step 2: Insert Chuck Brush into chuck mechanism.**



**Step 3: Turn chuck brush in a Circular Motion, Back & Forth for 5 seconds.**

**Step 4: Lubricate chuck mechanism with 1 Drop of Lubricant, Repeating Step 1.**

**Step 5: Run Handpiece (with Bur or Shipping Post) for 5 seconds to expel debris.**


**Chuck Brush Is Reusable!**  
**Please Ensure Suitable Cleaning of Brush after Every Use.**

[www.mcshandpiece.com/accessories](http://www.mcshandpiece.com/accessories)

**Lubricate & Clean the Handpiece.**

**Step 1:** Remove handpiece from unit tubing, clean with **Alcohol**.

**Step 2: Free Spin Your Turbine:** With **Bur Inserted** and **Coupling Attached**, place **Lubricant (2 - 3 Drops)** into the air drive port on the **coupling**, (**As Indicated in the Diagram Below**). At the same time, **Rotate the Bur** between your fingers, **Back & Forth**, for a few seconds.



**Step 3:** Re-Connect handpiece to unit tubing, and **Run for 30 Seconds** to eliminate excess oil from the turbine.

**Step 4:** After Autoclaving the handpiece, allow to cool gradually for **30 Minutes to Room Temperature**.

**Step 5: Repeat Steps 2 & 3.**

**We Recommend:**

For More Information about Maintenance: [www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

T: (604) 575-8544 - TF: (888) 575-8544 - TXT: (604) 373-6747 - [info@mcshandpiece.com](mailto:info@mcshandpiece.com) - [www.mcshandpiece.com](http://www.mcshandpiece.com)

**Office Hours:** Monday – Friday; 8am – 2pm (PST) - **Office Closures:** Please Visit Our Website for Dates!

Page 2 of 3

**Soaking your Handpieces in Lubricant Overnight; at Least Once a Month.**

*This will help Remove Stuck on Debris & Saturate Bearings with Lubricant, for longer life.*

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

---

## Helpful Links:

We have lots of extremely useful information on our website:

[www.mcshandpiece.com](http://www.mcshandpiece.com)

For an In-depth Look at Handpiece Maintenance:

[www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

Cause & Effect of Excessive Pressure Applied to the Bur:

[www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur](http://www.mcshandpiece.com/cause-effect-excessivepressureappliedtobur)

Understanding High Speed Bearings:

[www.mcshandpiece.com/understanding-high-speed-bearings](http://www.mcshandpiece.com/understanding-high-speed-bearings)

How to Free Spin Your Turbine:

[www.mcshandpiece.com/free-spin](http://www.mcshandpiece.com/free-spin)

Handpiece Lubricating Procedure:

[www.mcshandpiece.com/lubricating-procedure](http://www.mcshandpiece.com/lubricating-procedure)

Chuck Mechanism Cleaning Procedure:

[www.mcshandpiece.com/chuck-brush-procedure](http://www.mcshandpiece.com/chuck-brush-procedure)

Specialized Maintenance Procedures:

Highspeed's: [www.mcshandpiece.com/highspeed-maintenance](http://www.mcshandpiece.com/highspeed-maintenance)

Low Speed's: [www.mcshandpiece.com/low-speed-maintenance](http://www.mcshandpiece.com/low-speed-maintenance)

Electric Air Powered Highspeed's: [www.mcshandpiece.com/electric-airpowered-maintenance](http://www.mcshandpiece.com/electric-airpowered-maintenance)

Air Motors: [www.mcshandpiece.com/air-motor-maintenance](http://www.mcshandpiece.com/air-motor-maintenance)

Accessories & Other Purchases:

[www.mcshandpiece.com/products](http://www.mcshandpiece.com/products)

---

For More Information about Maintenance: [www.mcshandpiece.com/maintenance](http://www.mcshandpiece.com/maintenance)

T: (604) 575-8544 - TF: (888) 575-8544 - TXT: (604) 373-6747 - [info@mcshandpiece.com](mailto:info@mcshandpiece.com) - [www.mcshandpiece.com](http://www.mcshandpiece.com)

**Office Hours:** Monday – Friday; 8am – 2pm (PST) - **Office Closures:** Please Visit Our Website for Dates!

Page 3 of 3